



APAC Confirmation kits for urine & oral fluid laboratory testing

Purpose of usage

The APAC Confirmation kits are used for specimen collection, packaging and transportation of samples to the toxicology laboratory.

Chain of custody

The Tamperproof specimen bag and security seals, which cover the collection jar & transportation box, allows the client to monitor eventual tampering with the specimen, while it is shipped from the site of collection to the laboratory.

Scope of supply – urine confirmation kit

- Two specimen jars with screw lids
- Tamperproof specimen bag with adsorbent pad
- Four security seals; two to seal the specimen jars, one to seal the transport box and one to attach to the Chain of custody document. Each set of four seals uses a unique number for ease of tracking.
- Cardboard box to house specimen jars



Scope of supply – oral fluid confirmation kit

- Two specimen tubes and collection pads. The pads are used to collect oral fluid and will change color once sufficient fluid has been collected. The pad is inserted into the tube containing a saline solution, where after the tube is closed with a rubber cork.
- Tamperproof specimen bag with adsorbent pad
- Four security seals; two to seal the specimen tubes, one to seal the transport box and one to attach to the Chain of custody document. Each set of four seals uses a unique number for ease of tracking.
- Cardboard box to house specimen tubes



Keep spatulas in the mouth until colour has changed

You Configure - We Deliver

Call for a free sample

t: +61 2 9986 2252

e: info@apacsecurity.com

w: apacsecurity.com

